

SMART STRIP AFLA

**Lateral flow test kit for the detection
of total aflatoxins**
cod. MA810

SMART STRIP AFLA is a test kit for the quantitative and qualitative detection of total aflatoxins in grain commodities. The kit contains all necessary reagents and procedures to run 20 determinations.

Assay principle: Immunochromatographic competitive assay on strip. Digital analysis of results is mandatory

Analysable materials: corn

Sample preparation: grinding, methanol-water extraction, filtration, dilution with buffer

Qualitative analysis

Assay time: 5 minutes

Cut off: 10 ppb

Quantitative analysis

Assay time: 10 minutes

Measuring range: 2 – 75 ppb (10 – 375 ppb by sample dilution)

1 PROVIDED MATERIALS

Test strips: 2 aluminium bags containing 10 strips each

Dilution buffer: 1 plastic bottle containing 6 ml of ready to use buffer

Datamatrix: 1, batch specific

2 REQUIRED NOT PROVIDED MATERIALS

2.1 For sample preparation and assay implementation

- Distilled water and methanol; alternatively Mycotoxin Extraction Solution A, cod. ME070
- Grinder
- Balance
- For each sample, two plastic or glass container having at least 50 ml capacity (i.e. centrifuge tube, jar) and one further small vial for the final dilution (i.e. micro-centrifuge tube). One Whatman 1 filter.
- 50-200 μ l micropipette with proper tips
- 250 ml graduated cylinder, shaker, mini centrifuge - 5000 x g speed required- (optional)

2.2 For digital analysis

- Computer (Windows 7 32 or 64 bits, Windows 8 64 bits, Windows 8.1 64 bits, Windows 10 64 bits) with at least one USB port 2.0 or 3.0 (cod. AT600)
- LAB LFD reader (cod. LT900)

- Software Smart-Soft (cod. LT910)
- Template Smart-Strip- frame A (LT903)

3 PRECAUTIONS

- Store the kit at +2/+8°C, do not freeze.
- Take the strips from the bag and reseal immediately the unused cassettes together with provided desiccant.
- Don't use the kit after expiration.
- Don't use photocopies of the instruction booklet; use the kit insert only.
- For in vitro diagnostic use only.

4 SAMPLES PREPARATION

1. Mix carefully the sample to get it homogeneous.
2. Finely grind the sample.
3. Weight 50 g of ground sample and add 150 ml of 70% methanol in distilled water solution – or 150 ml of ready to use Mycotoxin Extraction Solution, cod. ME070. Alternatively: weigh 10 g of ground sample and add 30 ml of 70% methanol in distilled water solution.
4. Shake thoroughly 3 minutes. It is suggested to shake manually or by a mild shaker system.
5. Filter the sample on Whatman 1 and collect the filtrate. Alternatively: allow the sample to settle down then take 1 ml of the upper phase and centrifuge 1 minute at 5000 x g.
6. Dilute the extract 1:3 with the dilution buffer (100 µl of sample+ 200 µl of buffer). The sample is ready to be run in the strip and analysed in the range 2 – 75 ppb. *When calculating results, select “1:1” dilution factor in the software.*
7. For > 75 ppb samples, dilute the filtrate 1:5 with 70% methanol and then 1:3 in the dilution buffer. The measuring range is hence 10 – 375 ppb. *When calculating results, select “1:5” dilution factor in the software.*
8. Mix by pipetting up and down the liquids.

5 ASSAY PROCEDURE


5.1. Preliminary comments

- Bring all the kit components to room temperature prior opening the bag. Refer to point 3 “Precautions” relative to the proper use of strips. Return all reagents to +2/+8°C immediately after use.
- Get the LAB LFD reader ready.

- Don't change the assay protocol, don't modify the incubation time nor the sample preparation protocol.
- Use a single disposable tip for each sample to avoid any cross-contamination.
- The assay is not suitable for visual interpretation.

5.2 Assay procedure

1. In order to maintain sample traceability, label every lateral flow strip with sample code.
2. Take 100 µl of the diluted extract avoiding bubbles or foam and add it into the sample well of the cassette.
3. Incubate 5 or 10 minutes. The sample must be completely absorbed on the device.
4. When the incubation time expires, two coloured lines should be visible. One is labelled as T line (test line) and the other is C line (control line).
5. Read the developed strip within 1 minute. The colour intensity of C and T lines depends on the sample contamination. The C line is not supposed to disappear and should not be shaded. If C line is absent or smeared, consider the test as invalid. Do not put invalid tests into the reader.



6 RESULTS ACQUISITION

6.1 Qualitative analysis with LAB LFD Reader

1. Scan the Datamatrix and proceed with results acquisition.
2. When **5 minutes** incubation expires, put the developed cassette face down in one frame slot.
3. **For < 10 ppb samples the concentration result will be displayed as "<2 ppb". Higher contaminations than 10 ppb in the sample will be displayed as digits.** Such concentration does not correlate to the quantitative result and is not to be considered the sample concentration. **In order to determine the total Afla concentration, wait until 10 minutes incubation** and repeat the reading.

6.2 Quantitative analysis with LAB LFD reader

1. Scan the Datamatrix and proceed with results acquisition.
2. When **10 overall minutes** incubation expires, put the developed cassette face down in one frame slot.

WARNING: substitution will be possible just in case of rendered kit. The kit must be conserved in its integral version as indicated in this booklet.

7 LIABILITY

Those samples that are detected as positive have to be tested again with a confirmatory method. Tecna shall not be liable for any damage to the customer caused by the improper use of this kit and for any undertaken action as a consequence of results. Tecna shall not be liable for unsafe use of the kit out of the current European safety regulations.

The line for rapid screening of mycotoxins in cereals				
SMART STRIP AFLA	SMART STRIP AFLA B ₁	SMART STRIP DON	SMART STRIP FUMO	SMART STRIP ZON

DISTRIBUYE: