

Scope of delivery – ready-to-use

The BactoSonic® includes an ultrasonic bath BS 14, different sizes of implant containers, corresponding holders, other accessories and standard operating procedure of the sonication.



Endoprostheses not included.

The equipment is based on research and development cited in literature.

Akralab, S.L.
Pol. Industrial las Atalayas
Avda. de la Antigua Peseta - 77
Buzón 20212 - C.P: 03114 - Alicante
Tlf: 902 22 22 75 - Fax: 902 15 41 65
atencion.clientes@akralab.es
www.akralab.es

References

1. Achermann Y, Vogt M, Leunig M, Wust J, Trampuz A. Improved diagnosis of periprosthetic joint infection by multiplex PCR of sonication fluid from removed implants. J Clin Microbiol 2010; 48 (in press).
2. Bjerkan G, Witsø E, Bergh K. Sonication is superior to scraping for retrieval of bacteria in biofilm on titanium and steel surfaces in vitro. Acta Orthop 2009; 80:245-50.
3. Carmen JC, Roeder BL, Nelson JL, Ogilvie RL, Robison RA, Schaalje GB, Pitt WG. Treatment of biofilm infections on implants with low-frequency ultrasound and antibiotics. Am J Infect Control 2005; 33: 78-82.
4. Clauss M, Trampuz A, Borens O, Böhner M, Ilchmann T. Biofilm formation on bone grafts and bone graft substitutes. Comparison of different materials by a standard in vitro test and microcalorimetry. Acta Biomaterialia 2010; 6 (in press).
5. Del Pozo JL, Tran NV, Petty PM, Johnson CH, Walsh MF, Bite U, Clay RP, Mandrekar JN, Piper KE, Steckelberg JM, Patel R. Pilot study of association of bacteria on breast implants with capsular contracture. J Clin Microbiol 2009; 47: 1333-1337.
6. Kobayashi N, Bauer TW, Tuohy MJ, Fujishiro T, Procop GW. Brief ultrasonication improves detection of biofilm-formative bacteria around a metal implant. Clin Orthop Relat Res 2007;457: 210-3.
7. Monsen T, Lövgren E, Widerström M, Wallinder L. In vitro effect of ultrasound on bacteria and suggested protocol for sonication and diagnosis of prosthetic infections. J Clin Microbiol 2009; 47: 2496-501.
8. Piper KE, Jacobson MJ, Cofield RH, Sperling JW, Sanchez-Sotelo J, Osmon DR, Steckelberg JM, Mandrekar JN, Fernandez SM, Patel R. Microbiologic diagnosis of prosthetic shoulder infection using implant sonication. J Clin Microbiol 2009; 47: 1878-84.
9. Rieger UM, Pierer G, Lüscher NG, Trampuz A. Sonication of removed breast implants for improved detection of subclinical infection. Aesth Plast Surg 2009; 33: 404-408.
10. Rohacek M, Weissner M, Kobza R, Schoenenberger AW, Pfyffer G, Frei R, Erne P, Trampuz A. Bacterial colonization and infection of electrophysiologic cardiac devices detected with sonication and swab culture. Circulation 2010; 121 (in press).
11. Trampuz A, Osmon DR, Hanssen AD, Steckelberg JM, Patel R. Molecular and antibiofilm approaches to prosthetic joint infection. Clin Orthop Relat Res 2003; 69-88.
12. Trampuz A, Piper KE, Hanssen AD, Osmon DR, Cockerill FR, Steckelberg JM, Patel R. Sonication of explanted prosthetic components in bags for diagnosis of prosthetic joint infection is associated with a risk of contamination. J Clin Microbiol 2006; 44: 628-631.
13. Trampuz A., Piper KE, Jacobson MJ, Hanssen AD, Unni KK, Osmon DR, Mandrekar JN, Cockerill FR, Steckelberg JM, Greenleaf JF, Patel R. Sonication of removed hip and knee prostheses for diagnosis of infection. N Engl J Med 2007; 357: 654-663.
14. Trampuz A, Steinrücken J, Clauss M, Bizzini A, Furustrand U, Uckay I, Peter R, Bille J, Borens O. New methods for the diagnosis of implant-associated infections [article in French]. Rev Med Suisse 2010 (in press).
15. Zimmerli W, Trampuz A, Ochsner PE. Prosthetic-joint infections. N Engl J Med 2004; 351: 1645-1654.

BactoSonic®

Sonication

– optimized diagnostic method for implant-associated infections



Sz.: 6642e/2010-03

Implants improves the quality of life

Modern medicine has developed various **implants** to replace missing anatomical structure or biological function: joint prostheses, internal fixation devices, vascular prostheses, cardiac pacemakers and defibrillators, dental implants, neurosurgical shunts and breast implants.

New devices are improved and optimized with regard to biocompatibility and functionality.

Biofilm infections

– a challenge of modern medicine

With growing use of **implants**, modern medicine is facing an increasing risk of infections. Microorganisms on implant surface form **biofilms**, what makes them difficult to detect by conventional methods such as periprosthetic tissue cultures. For successful treatment of these infections **accurate microbiological diagnosis** is crucial. Such biofilms consist of an **extracellular matrix** of polymerized polysaccharide, in which bacteria are embedded (Fig. 1).




Figure 1 Formation of biofilm on implant surface

The microorganisms in biofilm are transformed in low metabolic, stationary growth state. Over weeks to months, depending on the type of microorganism, implant material and host, a **complex three-dimensional structure** develops, which consists of nutrition channels and rudimentary communication system resembling multicellular organism. Free-living (planktonic) bacteria are killed by antibiotics and the host defense system, while adherent (biofilm) bacteria can survive and persist in the extracellular matrix of the biofilm (Fig. 2).




Figure 2 Bacteria in biofilm resist antibiotics and defense system

Sonication – a new diagnostic method

After removal in the operating room, implants are placed in the **air-tight container** and transported to the microbiological laboratory. After addition of Ringer's solution, the implant is processed by vortexing (30 seconds) and sonication (1 minute) to **dislodge** (planktonize) microorganism into the surrounding fluid (sonicate). The sonication fluid is cultured on **aerobic and anaerobic agar plates** (Fig. 3) and inoculated in **broth media**.




Figure 3 Sonication removes more than 99.9 % of the biofilm bacteria from the surface

Principle of sonication

High acoustic intensity of conventional ultrasound baths kills microorganisms (especially gram-negative and anaerobic bacteria). Sonication in the specially designed ultrasound bath BactoSonic® uses **low frequency** and **low intensity ultrasound** at the threshold of microbubble formation (cavitation). Due micro-currents of sonication fluid, shear forces and oscillating cavitation bubbles biofilm is removed and the bacteria are disaggregated.




Figure 4 Sonication removes the detection of bacteria up to 10,000 times compared to periprosthetic tissue cultures

The resulting cavitation energy is reduced to the level, that no significant cell destruction occurs, enabling culture of viable microorganisms.

Advantages of sonication

• High accuracy

With low-intensity sonication, microorganisms are dislodged from the implant but not killed, enabling high sensitivity of conventional cultures (Fig. 4). Particularly **difficult to detect microorganisms** (including small-colony variants), **individual morphotypes** and **mixed infections** can be better detected. The sensitivity is particularly improved in patients receiving **previous antibiotics**, due to better survival of bacteria in biofilm. Ultrasound reaches through surrounding fluid the whole implant surface, which is associated with high specificity.

• Rapid result

Sonication increases microbial growth by inducing micro-currents in the sonication fluid, thereby shortening the microbial detection time.

• Quantitative biofilm assessment

Since bacteria survive, but not replicate in the sonication fluid, quantitative assessment of removed biofilm is possible. The microbial density is expressed as number of colony-forming units (CFU) per ml of sonication fluid.

• Additional investigations

The sonication fluid contains high density of bacteria, making it suitable for further **microbial** (e.g. PCR, MALDI-TOF, microcalorimetry) and **immunological analyses** (e.g. determination of biomarkers, gene expression).